

CE 0068

Siluk
MICRO FILLER NEEDLE

CE 0068



Visionmed Ltd.
B-1601 SK V1 GL Metrocity, 128 Beobwon-ro, Songpa-gu,
Seoul, 05854 Republic of Korea
Tel. +82-2-3482-1181 Fax. +82-2-3482-1182
E-mail. info@vmed.co.kr Website. www.vmed.co.kr

Premium
Micro
Filler Needle



Premium Micro Filler Needle

A new generation of cannulas
Designed and manufactured to
The highest levels of quality

“
Just feel the
difference
”



WHY Silk MICRO FILLER NEEDLE



Zoom

Precise Injection

Electro polished ultra thin wall cannula with
blunt tip enables as traumatic treatment
with all fillers

CE 0068

WHAT Silk MICRO FILLER NEEDLE



24Kit/Box
[with puncture needle]

- ▶ Ultra Thin Wall
- ▶ Excellent Results
- ▶ Less Trauma
- ▶ Minimal Downtime
- ▶ High quality Stainless Steel

8 different sizes [lengths, diameters and internal calibration] are
carefully adapted according to the treated areas, the direction
and the depth of injection

SPECIFICATION

Model No.	Cannula	Punching needle
FTMN18-70	18Gx70mm	18Gx25.4mm
FTMN22-50-21	22Gx50mm	21Gx25.4mm
FTMN22-70-21	22Gx70mm	21Gx25.4mm
FTMN23-50	23GX50mm	23GX25.4mm
FTMN25-40-23	25Gx40mm	23Gx25.4mm
FTMN25-50-23	25Gx50mm	23Gx25.4mm
FTMN27-25-26	27Gx25mm	26Gx12.7mm
FTMN27-40-26	27Gx40mm	26Gx12.7mm
FTMN27-50-26	27Gx50mm	26Gx12.7mm
FTMN30-25-27	30Gx25mm	27Gx12.7mm